

Product datasheet for **HP211152**

DNAJC6 Human qPCR Primer Pair (NM_014787)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9829
Forward Sequence:	CTGGACAGTCAGGAGTGGGAAGA
Reverse Sequence:	GGTGCTCCAAATGGGTCAAAGG
ACCN:	NM_014787 , NM_014787.1 , NM_014787.2 , NM_014787.3 , BC051722 , BC109279 , BC109280 , BP311649 , BQ423527 , BX647203 , NM_014787.4
UniProt ID:	O75061
Synonyms:	DJC6; PARK19
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)