

Product datasheet for **HP211163**

KIAA0528 (C2CD5) Human qPCR Primer Pair (NM_014802)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9847
Forward Sequence:	CTCAATCAGCGTTGGAACAGAGG
Reverse Sequence:	GTTTCTGGTTCATCAGGATTGTGG
ACCN:	<u>NM_014802</u> , <u>NM_014802.1</u> , <u>NM_014802.2</u> , <u>BC042498</u> , <u>BC042498.1</u> , <u>BC053885</u> , <u>BC117143</u> , <u>BC143859</u> , <u>BC143860</u> , <u>BC143878</u> , <u>BP418369</u> , <u>BX648560</u>
UniProt ID:	<u>Q86YS7</u>
Synonyms:	CDP138; KIAA0528
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)