

Product datasheet for **HP211168**

TMEM24 (C2CD2L) Human qPCR Primer Pair (NM_014807)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9854
Forward Sequence:	CTGTCTCGAAGACAGTTGTGCC
Reverse Sequence:	CTTCTTGGTAGGTGTGATGCTGG
ACCN:	<u>NM_014807</u> , <u>NM_014807.1</u> , <u>NM_014807.2</u> , <u>NM_014807.3</u> , <u>NM_014807.4</u> , <u>BC022219</u> , <u>BC022219.1</u> , <u>BC003598</u> , <u>BC010071</u> , <u>BC021254</u> , <u>BC033083</u> , <u>BC040131</u> , <u>BC052246</u> , <u>BC075832</u> , <u>BX327704</u> , <u>BX537572</u>
UniProt ID:	<u>O14523</u>
Synonyms:	DLNB23; TMEM24
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)