

## Product datasheet for **HP211174**

### PSMD6 Human qPCR Primer Pair (NM\_014814)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9861
Forward Sequence:	GAAGCGTTTGGATGAGGAGCTG
Reverse Sequence:	CTTGCGAAAGGCTGTCAGAGCT
ACCN:	<a href="#">NM_014814</a> , <a href="#">NM_014814.2</a> , <a href="#">BC012369</a> , <a href="#">BC012369.1</a> , <a href="#">BC000630</a> , <a href="#">BC000904</a> , <a href="#">BG723464</a> , <a href="#">BM681873</a> , <a href="#">BX647742</a> , <a href="#">NM_014814.3</a>
UniProt ID:	<a href="#">Q15008</a>
Synonyms:	p42A; p44S10; Rpn7; S10; SGA-113M
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)