

Product datasheet for **HP211195**

TECPR2 Human qPCR Primer Pair (NM_014844)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9895
Forward Sequence:	AGGAGATGACGACCATTGGTGG
Reverse Sequence:	GCGCAGCTTATCTGCCACCTTT
ACCN:	NM_014844 , NM_014844.1 , NM_014844.2 , NM_014844.3 , NM_014844.4 , BC030791 , BC136647 , BC142667 , BC142715 , BC151230 , BP365224 , BQ023338 , NM_014844.5
UniProt ID:	O15040
Synonyms:	KIAA0329; SPG49
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)