

Product datasheet for **HP211207**

RICH2 (ARHGAP44) Human qPCR Primer Pair (NM_014859)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9912
Forward Sequence:	AGCCAAGTTGGTGCTGGACATG
Reverse Sequence:	GCAAATCTCCACTCTGTTGGCAG
ACCN:	<u>NM_014859</u> , <u>NM_014859.1</u> , <u>NM_014859.2</u> , <u>NM_014859.3</u> , <u>NM_014859.5</u> , <u>BC022452</u> , <u>BC070074</u> , <u>BC094754</u> , <u>BC117412</u> , <u>BC117416</u> , <u>BC143852</u> , <u>BC143853</u> , <u>BC143855</u> , <u>BC143856</u> , <u>BM993036</u>
UniProt ID:	<u>Q17R89</u>
Synonyms:	NPC-A-10; RICH2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)