

Product datasheet for **HP211208**

STAF65 gamma (SUPT7L) Human qPCR Primer Pair (NM_014860)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9913
Forward Sequence:	CCTGATGACCTCCTGCCTTAG
Reverse Sequence:	GCTGAGTTCAGTCACAGGTCC
ACCN:	<u>NM_014860</u> , <u>NM_014860.1</u> , <u>NM_014860.2</u> , <u>BC074000</u> , <u>BC045638</u> , <u>BI755200</u> , <u>NM_014860.3</u>
UniProt ID:	<u>O94864</u>
Synonyms:	SPT7L; STAF65; STAF65(gamma); STAF65G; SUPT7H
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)