

## Product datasheet for **HP211243**

### Glutaminase (GLS) Human qPCR Primer Pair (NM\_014905)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2744
Forward Sequence:	CAGAAGGCACAGACATGGTTGG
Reverse Sequence:	GGCAGAAACCACCATTAGCCAG
ACCN:	<a href="#">NM_014905</a> , <a href="#">NM_014905.1</a> , <a href="#">NM_014905.2</a> , <a href="#">NM_014905.3</a> , <a href="#">NM_014905.4</a> , <a href="#">BC038507</a> , <a href="#">BC038507.1</a> , <a href="#">BM148085</a> , <a href="#">BX641148</a> , <a href="#">NM_014905.5</a>
UniProt ID:	<a href="#">O94925</a>
Synonyms:	AAD20; CASGID; DEE71; EIEE71; GAC; GAM; GDPAG; GLS1; KGA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)