

## Product datasheet for **HP211249**

### AAK1 Human qPCR Primer Pair (NM\_014911)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22848
Forward Sequence:	CAGCAGCAACTCTTCTCAAGC
Reverse Sequence:	GGTTGCTGGATGTACTGCCTGA
ACCN:	<u><a href="#">BC002695</a></u> , <u><a href="#">NM_014911</a></u> , <u><a href="#">NM_014911.1</a></u> , <u><a href="#">NM_014911.2</a></u> , <u><a href="#">NM_014911.3</a></u> , <u><a href="#">BC002695.2</a></u> , <u><a href="#">BC073953</a></u> , <u><a href="#">BC090950</a></u> , <u><a href="#">BC104842</a></u> , <u><a href="#">BC111965</a></u> , <u><a href="#">BC132937</a></u> , <u><a href="#">BC132939</a></u> , <u><a href="#">BC143710</a></u> , <u><a href="#">BX640963</a></u> , <u><a href="#">BX641060</a></u> , <u><a href="#">NM_014911.4</a></u>
UniProt ID:	<u><a href="#">Q2M218</a></u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)