

Product datasheet for **HP211253**

KPI2 (LMTK2) Human qPCR Primer Pair (NM_014916)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22853
Forward Sequence:	CGTAGCAAGAGTCATCGTGAAGG
Reverse Sequence:	GGAATCGCTTCTACGCACTGTC
ACCN:	<u>NM_014916</u> , <u>NM_014916.1</u> , <u>NM_014916.2</u> , <u>NM_014916.3</u> , <u>BC080576</u> , <u>BC131504</u> , <u>BC146773</u> , <u>BC167828</u> , <u>BE049530</u> , <u>BF001752</u> , <u>BM172505</u> , <u>NM_014916.4</u>
UniProt ID:	<u>Q8IWU2</u>
Synonyms:	AATYK2; BREK; cprk; hBREK; KPI-2; KPI2; LMR2; PPP1R100
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)