

Product datasheet for **HP211256**

ATG14L (ATG14) Human qPCR Primer Pair (NM_014924)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22863
Forward Sequence:	TGTACCTGGTCAGTCCAAGCTC
Reverse Sequence:	CAGGTCGGTTTCTTCATCGCTG
ACCN:	NM_014924 , NM_014924.1 , NM_014924.2 , NM_014924.3 , NM_014924.4 , BC109089 , BC028015 , BC109090 , BR000826 , NM_014924.5
UniProt ID:	Q6ZNE5
Synonyms:	ATG14L; BARKOR; KIAA0831
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)