

Product datasheet for **HP211259**

CNKSR2 Human qPCR Primer Pair (NM_014927)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22866
Forward Sequence:	CCACCAAAGAGACAGTTCTGCC
Reverse Sequence:	TATGCACTGGCTTGCCACCAT
ACCN:	<u>NM_014927</u> , <u>NM_014927.1</u> , <u>NM_014927.2</u> , <u>NM_014927.3</u> , <u>NM_014927.4</u> , <u>BC126121</u> , <u>BC092447</u> , <u>BC136289</u> , <u>BC143839</u> , <u>BM692904</u> , <u>BP200041</u> , <u>NM_014927.5</u>
UniProt ID:	<u>Q8WXI2</u>
Synonyms:	CNK2; KSR2; MAGUIN; MRXSHG
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)