

Product datasheet for **HP211282**

DIS3 Human qPCR Primer Pair (NM_014953)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22894
Forward Sequence:	GTAGACCCACCAGGATGTACTG
Reverse Sequence:	CACAGTTGTTCTTCTTGCTG
ACCN:	<u>NM_014953</u> , <u>NM_014953.1</u> , <u>NM_014953.2</u> , <u>NM_014953.3</u> , <u>NM_014953.4</u> , <u>BC056143</u> , <u>BC012468</u> , <u>BC038101</u> , <u>NM_014953.5</u>
UniProt ID:	<u>Q9Y2L1</u>
Synonyms:	2810028N01Rik; dis3p; EXOSC11; KIAA1008; RRP44
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)