

## Product datasheet for **HP211283**

### Rabphilin 3A (RPH3A) Human qPCR Primer Pair (NM\_014954)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22895
Forward Sequence:	GCCAACAGCTACGATTCGGATG
Reverse Sequence:	GGTGCAGCTTAACGTAGGGATC
ACCN:	<u><a href="#">NM_014954</a></u> , <u><a href="#">NM_014954.1</a></u> , <u><a href="#">NM_014954.2</a></u> , <u><a href="#">NM_014954.3</a></u> , <u><a href="#">BC017259</a></u> , <u><a href="#">BC017259.2</a></u> , <u><a href="#">BC002605</a></u> , <u><a href="#">BC009843</a></u> , <u><a href="#">BU738474</a></u> , <u><a href="#">BX647281</a></u>
UniProt ID:	<u><a href="#">Q9Y2J0</a></u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)