

## Product datasheet for **HP211287**

### CARD8 Human qPCR Primer Pair (NM\_014959)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22900
Forward Sequence:	GTGTGGGATACTGAGGTGAAGC
Reverse Sequence:	TGTCCTGGAGATCATCGAGCAC
ACCN:	<a href="#">NM_014959</a> , <a href="#">NM_014959.1</a> , <a href="#">NM_014959.2</a> , <a href="#">NM_014959.3</a> , <a href="#">BC056891</a> , <a href="#">BC045539</a> , <a href="#">BC034564</a> , <a href="#">BC041697</a> , <a href="#">BC046136</a> , <a href="#">BX641091</a> , <a href="#">BX648568</a> , <a href="#">NM_014959.5</a>
UniProt ID:	<a href="#">Q9Y2G2</a>
Synonyms:	CARDINAL; DACAR; DAKAR; NDPP; NDPP1; TUCAN
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)