

Product datasheet for **HP211304**

CEP152 Human qPCR Primer Pair (NM_014985)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22995
Forward Sequence:	GAGCCTTGGAAGAACTTCGTGG
Reverse Sequence:	TGCAGCTCGTTCCTTACTCTCC
ACCN:	<u>NM_014985</u> , <u>NM_014985.1</u> , <u>NM_014985.2</u> , <u>NM_014985.3</u> , <u>BC117182</u> , <u>BC029603</u> , <u>BC069186</u> , <u>BC107721</u> , <u>BG941609</u> , <u>BX648822</u> , <u>NM_014985.4</u>
UniProt ID:	<u>O94986</u>
Synonyms:	MCPH4; MCPH9; SCKL5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)