

Product datasheet for **HP211309**

DAAM1 Human qPCR Primer Pair (NM_014992)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23002
Forward Sequence:	ATTCCTCAGCCCACAAATGCC
Reverse Sequence:	GCAGAGAAGGTTCTTCCAGGTC
ACCN:	NM_014992 , NM_014992.1 , NM_014992.2 , BC024781 , BC038428 , BC064999 , BP199650 , BQ050640 , BX247986 , BX648383
UniProt ID:	Q9Y4D1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)