

Product datasheet for **HP211319**

SEMCAP3 (PDZRN3) Human qPCR Primer Pair (NM_015009)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23024
Forward Sequence:	TCTGGATGACCTGCACATGGAC
Reverse Sequence:	TGTCCTTCTCGTGCTGGTTGGA
ACCN:	BC014432 , NM_015009 , NM_015009.1 , NM_015009.2 , BC014432.1 , BC150623 , BC152417
UniProt ID:	Q9UPQ7
Synonyms:	LN3; SEMACAP3; SEMCAP3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)