

Product datasheet for **HP211328**

ZNF292 Human qPCR Primer Pair (NM_015021)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23036
Forward Sequence:	TCACGTCAGCTCCAACAAGGAG
Reverse Sequence:	CAAAGAAAGTTGACGACACCTCTC
ACCN:	<u>NM_015021</u> , <u>NM_015021.1</u> , <u>BC016386</u> , <u>BC040712</u> , <u>BC043168</u> , <u>BC063594</u> , <u>BC172360</u> , <u>BX537746</u> , <u>NM_015021.3</u>
UniProt ID:	<u>O60281</u>
Synonyms:	bA39312.3; MRD63; MRD64; Nbla00365; ZFP292; Zn-15; ZN-16
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)