

Product datasheet for **HP211340**

ZSWIM8 Human qPCR Primer Pair (NM_015037)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23053
Forward Sequence:	CAACTGGAACAGCGAGACACCT
Reverse Sequence:	GATGTTCTGCCTCGCTCACTGT
ACCN:	NM_015037 , NM_015037.1 , NM_015037.2 , NM_015037.3 , BC036835 , BC040726 , BC063694 , BC111006 , BC117205 , BC137249 , BC137250 , BC151206 , BC156550 , BM193415
UniProt ID:	A7E2V4
Synonyms:	KIAA0913
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)