

Product datasheet for **HP211348**

Senataxin (SETX) Human qPCR Primer Pair (NM_015046)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23064
Forward Sequence:	CAGTCCATCTCAAGACCTGTGC
Reverse Sequence:	GAGTTGAGCCATTCTGTGCCAC
ACCN:	NM_015046 , NM_015046.1 , NM_015046.2 , NM_015046.3 , NM_015046.4 , NM_015046.5 , BC032600 , BC032622 , BC078166 , BC106017 , BC137350 , BM718142 , BX537849 , BX538166
UniProt ID:	Q7Z333
Synonyms:	ALS4; AOA2; bA479K20.2; SCAN2; SCAR1; Sen1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)