

Product datasheet for **HP211373**

SARM (SARM1) Human qPCR Primer Pair (NM_015077)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23098
Forward Sequence:	CCAGATGTCTTCATCAGCTACCG
Reverse Sequence:	CACTCTGGATGAGTTTGCCTCG
ACCN:	<u>NM_015077</u> , <u>NM_015077.1</u> , <u>NM_015077.2</u> , <u>NM_015077.3</u> , <u>BC040429</u> , <u>BC122860</u> , <u>BG820141</u> , <u>NM_015077.4</u>
UniProt ID:	<u>Q6SZW1</u>
Synonyms:	hSARM1; HsTIR; MyD88-5; SAMD2; SARM
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)