

Product datasheet for **HP211381**

SPART Human qPCR Primer Pair (NM_015087)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23111
Forward Sequence:	AGGCAAATGTCTGACCTTCGGC
Reverse Sequence:	CACATCAGTGCCAGAGGCTTCT
ACCN:	<u>NM_015087</u> , <u>NM_015087.1</u> , <u>NM_015087.2</u> , <u>NM_015087.3</u> , <u>NM_015087.4</u> , <u>BC047083</u> , <u>BC026284</u> , <u>BQ719410</u> , <u>BU622545</u> , <u>BX648004</u> , <u>NM_015087.5</u>
UniProt ID:	<u>Q8N0X7</u>
Synonyms:	SPG20; TAHCCP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)