

Product datasheet for **HP211392**

ATG2A Human qPCR Primer Pair (NM_015104)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23130
Forward Sequence:	CCTCTGTGAGACCAAGGATGAG
Reverse Sequence:	CCAGACAGAAGTAGCCAAGTCC
ACCN:	NM_015104 , NM_015104.1 , NM_015104.2 , BC110650 , BC008593 , BC027481 , BC033850 , BC053596 , BC113090 , BC113091 , NM_015104.3
UniProt ID:	Q2TAZ0
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)