

Product datasheet for **HP211393**

RAD54L2 Human qPCR Primer Pair (NM_015106)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23132
Forward Sequence:	ACCACCTGGCTGTTCTTCTCT
Reverse Sequence:	TCCAGGAAGCATCAAACACCACC
ACCN:	<u>NM_015106</u> , <u>NM_015106.1</u> , <u>NM_015106.2</u> , <u>NM_015106.3</u> , <u>BC001474</u> , <u>BC024298</u> , <u>BQ014539</u> , <u>NM_015106.4</u>
UniProt ID:	<u>Q9Y4B4</u>
Synonyms:	ARIP4; HSPC325; SRISNF2L
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)