

Product datasheet for **HP211423**

RGL1 Human qPCR Primer Pair (NM_015149)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23179
Forward Sequence:	GCAAACCTGGACAGCAGTGTGA
Reverse Sequence:	AGTCCTGAAGGGCAGTGTCAAG
ACCN:	NM_015149 , NM_015149.1 , NM_015149.2 , NM_015149.3 , BC039250 , BC136591 , BF197850 , BP235344 , BX537416 , NM_015149.5
UniProt ID:	Q9NZL6
Synonyms:	RGL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)