

Product datasheet for **HP211452**

MAST4 Human qPCR Primer Pair (NM_015183)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	375449
Forward Sequence:	GATGGTGAGGCGGAGCAAGAAA
Reverse Sequence:	CAAGCAGGAGAACTTCGGCTG
ACCN:	<u>BC033215</u> , <u>NM_015183</u> , <u>NM_015183.1</u> , <u>NM_015183.2</u> , <u>BC033215.1</u> , <u>BC003646</u> , <u>BC018692</u> , <u>BF515864</u> , <u>BM967134</u> , <u>BX648796</u> , <u>BX648980</u> , <u>NM_015183.3</u>
UniProt ID:	<u>O15021</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)