

Product datasheet for **HP211464**

Cordon bleu homolog (COBL) Human qPCR Primer Pair (NM_015198)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23242
Forward Sequence:	CCAGTCCCAAAAGACACATCCC
Reverse Sequence:	TCCGTGTAGGACAGTTTCGCTG
ACCN:	NM_015198 , NM_015198.1 , NM_015198.2 , NM_015198.3 , NM_015198.4 , BC029275 , BC045771 , BC094695 , BC111496 , BC136441 , BC144099 , BC150263 , NM_015198.5
UniProt ID:	O75128
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)