

Product datasheet for **HP211488**

KIBRA (WWC1) Human qPCR Primer Pair (NM_015238)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23286
Forward Sequence:	GAACTGAGCCTTGGTAACAGCG
Reverse Sequence:	TGTCTCCAGCCACTGACTCGTC
ACCN:	NM_015238 , NM_015238.1 , NM_015238.2 , BC004394 , BC017746 , BC038463 , BC042553 , BC156527 , BC172519 , BX640827 , NM_015238.3
UniProt ID:	Q8IX03
Synonyms:	HBEBP3; HBEBP36; KIBRA; MEMRYQTL; PPP1R168
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)