

Product datasheet for **HP211503**

FAM21C (WASHC2C) Human qPCR Primer Pair (NM_015262)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	253725
Forward Sequence:	GAAGACCCAGAGAGTGCTCACTC
Reverse Sequence:	CATCGCTGAACAGCAAAGGTGG
ACCN:	NM_015262 , NM_015262.1 , NM_015262.2 , BC006456 , BC039561 , BC053887 , BC150611 , BM126842 , BQ229875 , NM_015262.3
UniProt ID:	Q9Y4E1
Synonyms:	FAM21A; FAM21C; VPEF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)