

Product datasheet for **HP211508**

DNAJC13 Human qPCR Primer Pair (NM_015268)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23317
Forward Sequence:	GGTTCCAATGTGCTTCCTGTTGC
Reverse Sequence:	GCTTCAGGTAGAAATGTGCCCAAG
ACCN:	NM_015268 , NM_015268.1 , NM_015268.2 , NM_015268.3 , BC151246 , BC000164 , BC009630 , BC033782 , BC040638 , BC043583 , BC065756 , BC105580 , NM_015268.4
UniProt ID:	O75165
Synonyms:	PARK21; RME8
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)