

Product datasheet for **HP211511**

RPGRIP1L Human qPCR Primer Pair (NM_015272)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23322
Forward Sequence:	GCCTCTGATTTTCGTTGGCACATG
Reverse Sequence:	GATCCACTTGGTGAAGGTAAGC
ACCN:	NM_015272 , NM_015272.1 , NM_015272.2 , NM_015272.3 , NM_015272.4 , BC017977 , BC136433 , BM683568 , BQ576116 , BX340768
UniProt ID:	Q68CZ1
Synonyms:	COACH3; CORS3; FTM; JBTS7; MKS5; NPHP8; PPP1R134
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)