

Product datasheet for **HP211516**

SASH1 Human qPCR Primer Pair (NM_015278)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23328
Forward Sequence:	GTAACAGCGACCAGTCAGGATC
Reverse Sequence:	GACTTGGCAGATAGGAGGCTAG
ACCN:	NM_015278 , NM_015278.1 , NM_015278.2 , NM_015278.3 , NM_015278.4 , BC028303 , BC063279 , BN000088 , BX537611 , NM_015278.5
UniProt ID:	O94885
Synonyms:	CAPOK; dj323M4.1; DUH1; SH3D6A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)