

## Product datasheet for **HP211523**

### PHF15 (JADE2) Human qPCR Primer Pair (NM\_015288)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23338
Forward Sequence:	CGAGAGAACGAAACACGCCATC
Reverse Sequence:	GCTGTTGAAGAAGGTGCCATCG
ACCN:	<a href="#">NM_015288</a> , <a href="#">NM_015288.1</a> , <a href="#">NM_015288.2</a> , <a href="#">NM_015288.3</a> , <a href="#">NM_015288.4</a> , <a href="#">NM_015288.5</a> , <a href="#">BC021962</a> , <a href="#">BC021962.2</a> , <a href="#">BC004292</a> , <a href="#">BC009307</a> , <a href="#">BC021169</a> , <a href="#">BN000288</a>
UniProt ID:	<a href="#">Q9NQC1</a>
Synonyms:	JADE-2; PHF15
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)