

Product datasheet for **HP211531**

KIAA0323 (KHNYN) Human qPCR Primer Pair (NM_015299)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23351
Forward Sequence:	TGTCATTGACGGCAGCAACGTG
Reverse Sequence:	CTTTGGCATCCTTACTGAAGCGC
ACCN:	NM_015299 , NM_015299.1 , NM_015299.2 , BC042554 , BC042554.1 , BC015780 , BM464048 , BX161465 , NM_015299.3
UniProt ID:	O15037
Synonyms:	KIAA0323
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)