

Product datasheet for **HP211597**

HECTD1 Human qPCR Primer Pair (NM_015382)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	25831
Forward Sequence:	GAGAAGGCACTGTCACAGGAGA
Reverse Sequence:	GATGCCACTGTATCAGGGTCGT
ACCN:	NM_015382 , NM_015382.1 , NM_015382.3 , BC006237 , BC011658 , BC016947 , BC063686 , BC072441 , BC083494 , BC172391 , BM977042
UniProt ID:	Q9ULT8
Synonyms:	EULIR
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)