

Product datasheet for **HP211601**

YIPF3 Human qPCR Primer Pair (NM_015388)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	25844
Forward Sequence:	CTTTGATGTGGAGCCTGCTCAG
Reverse Sequence:	GGATAGCAACCAGAGTGAAGACC
ACCN:	NM_015388 , NM_015388.1 , NM_015388.2 , NM_015388.3 , BC019297 , BC019297.1 , BE562740 , BM836904 , NM_015388.4
UniProt ID:	Q9GZM5
Synonyms:	C6orf109; dJ337H4.3; FinGER3; KLIP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)