

## Product datasheet for **HP211617**

### AASDHPPT Human qPCR Primer Pair (NM\_015423)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	60496
Forward Sequence:	GCAGTTGCTCTTAGGAAACCCG
Reverse Sequence:	AGGTGTCATGGGAACGGCAGAT
ACCN:	<a href="#">NM_015423</a> , <a href="#">NM_015423.2</a> , <a href="#">BC016728</a> , <a href="#">BC016728.1</a> , <a href="#">BC015470</a> , <a href="#">BX537665</a> , <a href="#">NM_015423.3</a>
UniProt ID:	<a href="#">Q9NRN7</a>
Synonyms:	AASD-PPT; ACPS; CGI-80; LYS2; LYS5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)