

Product datasheet for **HP211641**

ATL3 Human qPCR Primer Pair (NM_015459)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	25923
Forward Sequence:	TGGACTCCAAGGAGGAATGGCA
Reverse Sequence:	GAGAAAGCAGGTGACATCGGAG
ACCN:	NM_015459 , NM_015459.1 , NM_015459.2 , NM_015459.3 , NM_015459.4 , BC077727 , BC077727.1 , BE710391 , BI829120 , BX642720 , BX648468 , NM_015459.5
UniProt ID:	Q6DD88
Synonyms:	HSN1F
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)