

## Product datasheet for **HP211647**

### HDPTP (PTPN23) Human qPCR Primer Pair (NM\_015466)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	25930
Forward Sequence:	ACATCACTGCCTCGCTGGTCAC
Reverse Sequence:	GTTCTGCTCCAGGTACACCTTC
ACCN:	<a href="#">BC004881</a> , <a href="#">NM_015466</a> , <a href="#">NM_015466.1</a> , <a href="#">NM_015466.2</a> , <a href="#">NM_015466.3</a> , <a href="#">BC089042</a> , <a href="#">BC089042.1</a> , <a href="#">BC020380</a> , <a href="#">BC027711</a> , <a href="#">BT009758</a> , <a href="#">BU173747</a> , <a href="#">NM_015466.4</a>
UniProt ID:	<a href="#">Q9H3S7</a>
Synonyms:	HD-PTP; HDPTP; NEDBASS; PTP-TD14
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)