

Product datasheet for **HP211650**

TAZ (WWTR1) Human qPCR Primer Pair (NM_015472)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	25937
Forward Sequence:	GAGGACTTCCTCAGCAATGTGG
Reverse Sequence:	CGTTTGTTCTGGAAGACAGTCA
ACCN:	NM_015472 , NM_015472.1 , NM_015472.2 , NM_015472.3 , NM_015472.4 , NM_015472.5 , BC014052 , BC014052.2 , BQ003163 , NM_015472.6
UniProt ID:	Q9GZV5
Synonyms:	TAZ
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)