

Product datasheet for **HP211655**

mSin3A (SIN3A) Human qPCR Primer Pair (NM_015477)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	25942
Forward Sequence:	CAGAATGACACCAAGGTCCTGAG
Reverse Sequence:	CATACGCAAGTGAGAGGTGTGG
ACCN:	<u>NM_015477</u> , <u>NM_015477.1</u> , <u>NM_015477.2</u> , <u>BC137098</u> , <u>BC018973</u> , <u>BC030563</u> , <u>BC041047</u> , <u>BC052275</u> , <u>BC066364</u> , <u>BC137099</u> , <u>BM676304</u> , <u>BX649162</u> , <u>NM_015477.3</u>
UniProt ID:	<u>Q96ST3</u>
Synonyms:	WITKOS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)