

## Product datasheet for **HP211673**

### NECAP1 Human qPCR Primer Pair (NM\_015509)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	25977
Forward Sequence:	ATTGGACTGGTCGCCTCCGAAT
Reverse Sequence:	TAGCGGCTAGAATCTGTACCG
ACCN:	<a href="#">BD127364</a> , <a href="#">NM_015509</a> , <a href="#">NM_015509.1</a> , <a href="#">NM_015509.2</a> , <a href="#">NM_015509.3</a> , <a href="#">BC110876</a> , <a href="#">BC002888</a> , <a href="#">BC067367</a> , <a href="#">BC084551</a> , <a href="#">NM_015509.4</a>
UniProt ID:	<a href="#">Q8NC96</a>
Synonyms:	DEE21; EIEE21
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)