

## Product datasheet for **HP211681**

### IBTK Human qPCR Primer Pair (NM\_015525)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	25998
Forward Sequence:	GAAGGACAATCATTGGCGTTGCA
Reverse Sequence:	AGAGACCTGACGAGGAGCAGTT
ACCN:	<u><a href="#">NM_015525</a></u> , <u><a href="#">NM_015525.1</a></u> , <u><a href="#">NM_015525.2</a></u> , <u><a href="#">NM_015525.3</a></u> , <u><a href="#">BC113696</a></u> , <u><a href="#">BC012551</a></u> , <u><a href="#">BC027483</a></u> , <u><a href="#">BC027490</a></u> , <u><a href="#">BC035860</a></u> , <u><a href="#">BC038244</a></u> , <u><a href="#">BC042171</a></u> , <u><a href="#">BC092410</a></u> , <u><a href="#">BC112929</a></u> , <u><a href="#">BC113698</a></u> , <u><a href="#">BC143736</a></u> , <u><a href="#">BC143745</a></u> , <u><a href="#">BC143748</a></u> , <u><a href="#">BM670666</a></u> , <u><a href="#">NM_015525.4</a></u>
UniProt ID:	<u><a href="#">Q9P2D0</a></u>
Synonyms:	BTBD26; BTKI
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)