

Product datasheet for **HP211712**

AUTS2 Human qPCR Primer Pair (NM_015570)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26053
Forward Sequence:	CCTCCTCATCACAGCAACTTC
Reverse Sequence:	GAAGGCATTGCCACCAACTGCT
ACCN:	<u>NM_015570</u> , <u>NM_015570.1</u> , <u>NM_015570.2</u> , <u>NM_015570.3</u> , <u>BC064693</u> , <u>BC011643</u> , <u>BG208279</u> , <u>BU608023</u> , <u>NM_015570.4</u>
UniProt ID:	<u>Q8WXX7</u>
Synonyms:	FBRSL2; MRD26
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)