

## Product datasheet for **HP211719**

### GPATCH4 Human qPCR Primer Pair (NM\_015590)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54865
Forward Sequence:	TGAGATGCTGCTCCAAGCCTGT
Reverse Sequence:	GCCTTTGAGACGAGCCAGGAAG
ACCN:	<u><a href="#">NM_015590</a></u> , <u><a href="#">NM_015590.2</a></u> , <u><a href="#">NM_015590.3</a></u> , <u><a href="#">BC032436</a></u> , <u><a href="#">BC040147</a></u> , <u><a href="#">BC056904</a></u> , <u><a href="#">BC152774</a></u> , <u><a href="#">BE326312</a></u> , <u><a href="#">BF669837</a></u> , <u><a href="#">BM997996</a></u> , <u><a href="#">BQ214732</a></u>
Synonyms:	GPATC4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)