

## Product datasheet for **HP211723**

### GPSM1 Human qPCR Primer Pair (NM\_015597)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26086
Forward Sequence:	CCTGTTGACCAAGTTCCAGAGC
Reverse Sequence:	GGTCGAAGAATTCCTCGGTCTG
ACCN:	<u><a href="#">BC017353</a></u> , <u><a href="#">NM_015597</a></u> , <u><a href="#">NM_015597.1</a></u> , <u><a href="#">NM_015597.2</a></u> , <u><a href="#">NM_015597.3</a></u> , <u><a href="#">NM_015597.4</a></u> , <u><a href="#">NM_015597.5</a></u> , <u><a href="#">BC017353.1</a></u> , <u><a href="#">BC009943</a></u> , <u><a href="#">BC009979</a></u> , <u><a href="#">BC017257</a></u> , <u><a href="#">BC048343</a></u> , <u><a href="#">BI766024</a></u> , <u><a href="#">BQ214873</a></u> , <u><a href="#">BU631041</a></u>
UniProt ID:	<u><a href="#">Q86YR5</a></u>
Synonyms:	AGS3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)