

Product datasheet for **HP211727**

WDR21A (DCAF4) Human qPCR Primer Pair (NM_015604)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26094
Forward Sequence:	ACCACGAAGTGCCTAAGGCAGT
Reverse Sequence:	AGGCTCCAGATTCTCGTGTAGC
ACCN:	<u>NM_015604</u> , <u>NM_015604.1</u> , <u>NM_015604.2</u> , <u>NM_015604.3</u> , <u>BC018979</u> , <u>BC035901</u> , <u>BC035901.1</u> , <u>BX161489</u> , <u>BX247962</u>
UniProt ID:	<u>Q8WW16</u>
Synonyms:	WDR21; WDR21A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)