

## Product datasheet for **HP211729**

### **C10orf137 (EDRF1) Human qPCR Primer Pair (NM\_015608)**

#### **Product data:**

|                          |   |
|--------------------------|---|
| <b>Product Type:</b>     | qPCR Primer Pairs   |
| <b>Gene ID:</b>          | 26098   |
| <b>Forward Sequence:</b> | TTGTCCGATGCTGCCATGAGTC  |
| <b>Reverse Sequence:</b> | CGGAAAGCATATTGGTGAGAGG  |
| <b>ACCN:</b>             | <u><a href="#">BC026172</a></u> , <u><a href="#">BC029650</a></u> , <u><a href="#">NM_015608</a></u> , <u><a href="#">NM_015608.1</a></u> , <u><a href="#">NM_015608.2</a></u> , <u><a href="#">BC026172.1</a></u> , <u><a href="#">BC016903</a></u> ,<br><u><a href="#">BC028715</a></u> , <u><a href="#">BC105929</a></u> , <u><a href="#">BC107479</a></u> , <u><a href="#">BC107896</a></u> |
| <b>UniProt ID:</b>       | <u><a href="#">Q3B7T1</a></u>   |
| <b>Synonyms:</b>         | C10orf137   |
| <b>Components:</b>       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)  |
| <b>Quality Control:</b>  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.  |
| <b>Storage:</b>          | Store at -20°C.   |
| <b>Stability:</b>        | The primer mix is stable for one year from date of shipping.  |



[View online »](#)