

Product datasheet for **HP211737**

PRPF31 Human qPCR Primer Pair (NM_015629)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26121
Forward Sequence:	GCTTCTCGTCTACCTCAGTGCT
Reverse Sequence:	TGTGCTCTCGTGAAACTGTCC
ACCN:	NM_015629 , NM_015629.2 , NM_015629.3 , BC117389 , NM_015629.4
UniProt ID:	Q8WWY3
Synonyms:	NY-BR-99; PRP31; RP11; SNRNP61
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)